



# BrailleNote Touch Command Summary

The Command summary chapter provides a list of commands used in all applications.

## General

Home button (Main Menu)	<i>SPACE</i> with <i>DOTS</i> 1 2 3 4 5 6
Back button	<i>SPACE</i> with <i>E</i>
Enter	<i>DOT 8</i>
TTS Volume up	<i>ENTER</i> with <i>DOT 4</i>
TTS Volume down	<i>ENTER</i> with <i>DOT 1</i>
Notification volume up	<i>ENTER</i> with <i>DOT 5</i>
Notification volume down	<i>ENTER</i> with <i>DOT 2</i>
Increase speech rate	<i>ENTER</i> with <i>DOT 6</i>
Decrease speech rate	<i>ENTER</i> with <i>DOT 3</i>
Open KeySoft Options	<i>SPACE</i> with <i>O</i>
Open screen describer	<i>BACKSPACE</i> with <i>H</i>
Stop speech	<i>BACKSPACE</i> with <i>ENTER</i>
Repeat	<i>SPACE</i> with <i>R</i>
Activate notification shade	<i>ENTER</i> with <i>N</i>
Jump to any item in a list or on a home screen	Type the first letter of the desired item or app
Running apps	<i>SPACE</i> with dots 2-3-5 or press and hold the <b>Contextual menu button</b>
Clear app from app switcher	<i>BACKSPACE</i> with <i>C</i>
Contextual help	<i>SPACE</i> with <i>H</i>
Contextual menu	<i>SPACE</i> with <i>M</i> or <b>Contextual menu button</b>
Quick settings	<i>ENTER</i> with <i>Q</i>
System information	<i>SPACE</i> with <i>I</i>
Toggle TouchBraille On/off	<i>Previous</i> + <i>Next</i> thumb keys
Toggle language profile	<i>ENTER</i> with <i>L</i>
Radio status (Wi-Fi and Bluetooth)	<i>ENTER</i> with <i>DOTS</i> 3 4

Turn screen reader and accessibility on / off

Triple press the **Home button**

Unlock screen

**SPACE** with **U**

Toggle speech on/off

**SPACE** with **Previous** thumb key

Toggle visual display on/off

**SPACE** with **Next** thumb key

Eject media

**ENTER** with **E**

Toggle braille grade

**BACKSPACE** with **G**

## Navigation

Activate the selected item

**ENTER** or cursor routing button above the word on the display

Move to the previous item

**SPACE** with **DOT 1** or **Previous** thumb key

Move to the next item

**SPACE** with **DOT 4** or **Next** thumb key

Braille panning left and right

**Left** or **Right** thumb key

Move to first element on visual screen

**ENTER** with **DOTS 1 3**

Move to last element on visual screen

**ENTER** with **DOTS 4 6**

Scroll up in a list

**ENTER** with **DOTS 1 2 6**

Scroll down in a list

**ENTER** with **DOTS 3 4 5**

Toggle navigation levels

**SPACE** with **T** repeatedly

Move back using selected navigation element

**SPACE** with **DOTS 1 3**

Move forward using selected navigation element

**SPACE** with **DOTS 4 6**

Long press

**BACKSPACE** with **L**

Time

**ENTER** with **T**

Date

**ENTER** with **D**

Battery power

**ENTER** with **P**

# Reading

Previous character	<b>SPACE</b> with <b>DOT 3</b>
Current character	<b>SPACE</b> with <b>DOTS 3 6</b>
Next character	<b>SPACE</b> with <b>DOT 6</b>
Previous word	<b>SPACE</b> with <b>DOT 2</b>
Current word	<b>SPACE</b> with <b>DOTS 2 5</b>
Next word	<b>SPACE</b> with <b>DOT 5</b>
Previous line	<b>SPACE</b> with <b>DOT 1</b>
Current line	<b>SPACE</b> with <b>DOTS 1 4</b>
Next line	<b>SPACE</b> with <b>DOT 4</b>
Previous paragraph	<b>SPACE</b> with <b>DOTS 2 3</b>
Current paragraph	<b>SPACE</b> with <b>DOTS 2 3 5 6</b>
Next paragraph	<b>SPACE</b> with <b>DOTS 5 6</b>
Jump to the beginning of a line or sentence (Home)	<b>ENTER</b> with <b>DOTS 1 3</b>
Jump to the end of a line or sentence (End)	<b>ENTER</b> with <b>DOTS 4 6</b>
Jump to the top of an area	<b>SPACE</b> with <b>DOTS 1 2 3</b>
Jump to the bottom of an area	<b>SPACE</b> with <b>DOTS 4 5 6</b>
Read all	<b>SPACE</b> with <b>G</b>
Stop reading	<b>BACKSPACE</b> with <b>ENTER</b>

# Editing text

Activate edit mode	<b>ENTER</b> or cursor routing button
Leave edit mode	<b>SPACE</b> with <b>E</b> or Back
Move to next edit box while editing	<b>ENTER</b> or routing key
Move to next edit box without editing	<b>SPACE</b> with <b>DOT 4</b> or <b>Next</b> thumb key
Move to previous edit box without editing	<b>SPACE</b> with <b>DOT 1</b> or <b>Previous</b> thumb key
Move insertion point to start of text field document	<b>SPACE</b> with <b>DOTS 1 2 3</b> (jumps to the top of the edit text field)
Move insertion point to end of text field document	<b>SPACE</b> with <b>DOTS 4 5 6</b> (jumps to the end of the edit text field)
Delete previous character	<b>BACKSPACE</b> ( <b>DOT 7</b> )
Delete previous word	<b>BACKSPACE</b> with <b>DOT 2</b>
Delete current word	<b>BACKSPACE</b> with <b>DOTS 2 5</b>
Delete from cursor to the end of the line	<b>BACKSPACE</b> with <b>C</b>
Delete from cursor to the end of the document	<b>BACKSPACE</b> with <b>DOTS 4 5 6</b>
Start and stop selection	<b>ENTER</b> with <b>S</b>
Select all	<b>ENTER</b> with <b>DOTS 1 2 3 4 5 6</b>
Cut	<b>BACKSPACE</b> with <b>X</b>
Copy	<b>BACKSPACE</b> with <b>Y</b>
Paste	<b>BACKSPACE</b> with <b>V</b>

# KeyWord

Save	<b>SPACE</b> with <b>S</b>
Save as	<b>BACKSPACE</b> with <b>S</b>
File info	<b>SPACE</b> with <b>I</b>
Preview	<b>ENTER</b> with <b>V</b>
Spell check	<b>SPACE</b> with <b>DOTS 1 6</b>
Find	<b>SPACE</b> with <b>F</b>
Find previous	<b>SPACE</b> with <b>P</b>
Find next	<b>SPACE</b> with <b>N</b>
Toggle reading mode	<b>SPACE</b> with <b>X</b>
Tab	<b>BACKSPACE</b> with <b>T</b>
Align left	<b>BACKSPACE</b> with <b>L</b>
Align center	<b>ENTER</b> with <b>C</b>
Align right	<b>ENTER</b> with <b>R</b>
Align justify	<b>ENTER</b> with <b>J</b>
Bold	<b>ENTER</b> with <b>B</b>
Italic	<b>ENTER</b> with <b>I</b>
Underline	<b>ENTER</b> with <b>U</b>
Formatting details	<b>SPACE</b> with <b>DOTS 3 5</b>
Unicode Map	<b>BACKSPACE</b> with <b>DOTS 3 5</b>

# KeyFiles

Mark/Unmark

***BACKSPACE*** with ***L***

Mark/Unmark all

***ENTER*** with ***DOTS 1 2 3 4 5 6***

Show more info

***SPACE*** with ***I***

Rename

***BACKSPACE*** with ***R***

Delete

***BACKSPACE*** with ***DOTS 2 3 5 6***

Copy

***BACKSPACE*** with ***Y***

Cut

***BACKSPACE*** with ***X***

Paste

***BACKSPACE*** with ***V***

New folder

***SPACE*** with ***N***

Search

***SPACE*** with ***F***

Sort

***SPACE*** with ***V***

Where am I

***SPACE*** with ***Dots 1 5 6***

Select drive

***SPACE*** with ***D***

# KeyMail

Reply	<b>BACKSPACE</b> with <b>R</b>
Reply all	<b>BACKSPACE</b> with <b>A</b>
Forward	<b>BACKSPACE</b> with <b>F</b>
New message	<b>SPACE</b> with <b>N</b>
Mark a message in the message list	<b>BACKSPACE</b> With <b>L</b>
Send email	<b>BACKSPACE</b> with <b>S</b>
Refresh to check for new mail	<b>ENTER</b> with <b>R</b>
Delete email	<b>BACKSPACE</b> with <b>DOTS 2 3 5 6</b>
Save message as draft	<b>SPACE</b> with <b>S</b>
Search	<b>SPACE</b> with <b>F</b>
Change email folder	<b>ENTER</b> with <b>B</b>

# VictorReader

Bookshelf menu	<b><i>SPACE</i></b> with <b><i>B</i></b>
Go to Option menu	<b><i>ENTER</i></b> with <b><i>G</i></b>
Bookmark menu	<b><i>ENTER</i></b> with <b><i>M</i></b>
Navigation Level menu	<b><i>SPACE</i></b> with <b><i>T</i></b>
Move to previous element	<b><i>SPACE</i></b> with <b><i>DOTS 1 3</i></b>
Move to next element	<b><i>SPACE</i></b> with <b><i>DOTS 4 6</i></b>
Move to previous element while playing	<b><i>Previous</i></b> thumb key
Move to next element while playing	<b><i>Next</i></b> thumb key
Confirm	<b><i>Enter</i></b>
Start or Stop playback	<b><i>SPACE</i></b> with <b><i>G</i></b>
Stop playback	<b><i>ENTER</i></b> with <b><i>BACKSPACE</i></b>
Decrease speed	<b><i>ENTER</i></b> with <b><i>DOT 3</i></b>
Increase speed	<b><i>ENTER</i></b> with <b><i>DOT 6</i></b>
Toggle between playback modes	<b><i>ENTER</i></b> with <b><i>V</i></b>
Undo last navigation	<b><i>BACKSPACE</i></b> with <b><i>Z</i></b>
Go to beginning of book	<b><i>SPACE</i></b> with <b><i>DOTS 1 2 3</i></b>
Go to end of book	<b><i>SPACE</i></b> with <b><i>DOTS 4 5 6</i></b>
Sleep	<b><i>SPACE</i></b> with <b><i>S</i></b>
Rewind	Left press and hold while playing
Fast Forward	Right press and hold while playing



# KeyWeb

Open Web Page	<b>ENTER</b> with <b>O</b>
Back	<b>SPACE</b> with <b>E</b>
Forward	<b>ENTER</b> with <b>F</b>
Home	<b>ENTER</b> with <b>H</b>
Reload	<b>ENTER</b> with <b>R</b>
Bookmark Current Page	<b>ENTER</b> with <b>M</b>
Bookmarks	<b>ENTER</b> with lower <b>F</b>
History	<b>ENTER</b> with lower <b>H</b>
Reader Mode	<b>BACKSPACE</b> with <b>R</b>
Downloads	<b>SPACE</b> with <b>D</b>
New Tab	<b>BACKSPACE</b> with <b>N</b>
Tab List	<b>BACKSPACE</b> with <b>T</b>
Find	<b>SPACE</b> with <b>F</b>
Find next	<b>SPACE</b> with <b>N</b>
Find Previous	<b>SPACE</b> with <b>P</b>
Cycle Html Navigation Levels	<b>SPACE</b> with <b>T</b>
Next Html Element	<b>SPACE</b> with <b>DOTS 4 6</b>
Previous Html Element	<b>SPACE</b> with <b>DOTS 1 3</b>
Top of page	<b>SPACE</b> with <b>DOTS 1 2 3</b>
Bottom of page	<b>SPACE</b> with <b>DOTS 4 5 6</b>
Next block of text	<b>ENTER</b> with <b>DOTS 5 6</b>
Previous block of text	<b>ENTER</b> with <b>DOTS 2 3</b>
KeyWeb menu	<b>SPACE</b> with <b>W</b>

## KeyList

New contact	<b>SPACE</b> with <b>N</b>
Find contact	<b>SPACE</b> with <b>F</b>
Edit contact	<b>BACKSPACE</b> with <b>E</b>
Delete contact	<b>BACKSPACE</b> with <b>DOTS 2 3 5 6</b>
Save contact	<b>SPACE</b> with <b>S</b>
Refresh contacts	<b>ENTER</b> with <b>R</b>

## KeyPlan

New appointment	<b>SPACE</b> with <b>N</b>
Edit appointment	<b>BACKSPACE</b> with <b>E</b>
Delete appointment	<b>BACKSPACE</b> with <b>DOTS 2 3 5 6</b>
Go to date	<b>ENTER</b> with <b>G</b>
Today	<b>SPACE</b> with <b>T</b>
Where am I	<b>SPACE</b> with <b>DOTS 1 5 6</b>
Next day	<b>SPACE</b> with <b>DOT 6</b>
Previous day	<b>SPACE</b> with <b>DOT 3</b>
Next week	<b>SPACE</b> with <b>DOT 5</b>
Previous week	<b>SPACE</b> with <b>DOT 2</b>
Next month	<b>SPACE</b> with <b>DOT 4</b>
Previous month	<b>SPACE</b> with <b>DOT 1</b>
Next year	<b>SPACE</b> with <b>DOTS 5 6</b>
Previous year	<b>SPACE</b> with <b>DOTS 2 3</b>
Refresh appointments	<b>ENTER</b> with <b>R</b>

## KeyCalc US

Plus	<i><b>DOTS 3 4 6</b></i>
Minus	<i><b>DOTS 3-6</b></i>
Multiplication	<i><b>DOTS 1-6</b></i>
Divided by	<i><b>DOTS 3-4</b></i>
Equals	<i><b>ENTER</b></i>
Clear	<i><b>SPACE with DOTS 3-5-6</b></i>
Decimal point	<i><b>DOTS 4-6</b></i>
Percent	<i><b>DOTS 1-4-6</b></i>
Square root	<i><b>SPACE with DOTS 3-4-5</b></i>
Factorial	<i><b>SPACE with DOTS 2-3-4-6</b></i>
Power	<i><b>SPACE with DOTS 4-5</b></i>
Pi	<i><b>SPACE with Y</b></i>

## KeyCalc UK

Plus	<i><b>DOTS 2 3 5</b></i>
Minus	<i><b>DOTS 3 6</b></i>
Multiplication	<i><b>DOTS 3 5</b></i>
Divided by	<i><b>DOTS 1 2 3 6 8</b></i>
Equals	<i><b>ENTER</b></i>
Clear	<i><b>SPACE with DOTS 3 5 6</b></i>
Decimal point	<i><b>DOTS 2 5 6</b></i>
Percent	<i><b>DOTS 4 6</b></i>
Square root	<i><b>Dots 1 2 3 5 6 8</b></i>
Factorial	<i><b>SPACE with DOTS 2 3 4 6</b></i>
Power	<i><b>SPACE with DOTS 4 5</b></i>
Pi	<i><b>SPACE with Y</b></i>